

Relaxation values for BioPAL products

R1 and R2 values were determined using a Bruker Minispec mq20 (field 0.47T; temperature 37°C; solvent; water).

Product Code	Lot Number	R1, (mM ⁻¹ sec ⁻¹)	R2, (mM ⁻¹ sec ⁻¹)
--------------	------------	---	---

Molday ION			
CL-30Q02-2	09091101	36.4	70.6
CL-30Q02-2	10020401	--	62.2
CL-30Q02-7	11072201	32.4	62.7
CL-30Q02-CA	09112001	34.8	62.1
CL-30Q02-6A	10120601	32.1	67.1
CL-50Q02-6A-50	10062301	30.4	75.8

FeRex			
MR-7200	0509	25.2	160.6

Gd Based Products			
CL-20P02-2	09072501	6.8	7.4
CL-50P02-6	0505	4.8	5.2
P-00P01	07012401	9.1	10.0