



80 Webster Street
Worcester, MA 01603 (USA)
TEL: 508-770-1190

Web Site: BioPAL.com
e-mail: cs@BioPAL.com

Product Data Sheet

Product Name:	Molday ION Coumarin
Catalog Number:	CL-50Q02-6A-52
Description:	Stable aqueous suspension of iron-based superparamagnetic nanoparticles (USPIO)
Application:	Used to label cells, providing MRI contrast, magnetic targeting, capture and retention.
Packaging:	2 ml sealed serum bottle
Volume:	2.0 ml
Weight of Iron:	2 mg Fe/ml
Mean Particle Size:	35 nm
Core:	Magnetite
Mean Core Size:	8 nm
Coating:	Dextran
Number of Particles:	~6.4E16 per g
Density:	~1.25 g/cm ³
Zeta Potential:	~ +31 mV
R1 (mM ⁻¹ sec ⁻¹)	30.4
R2 (mM ⁻¹ sec ⁻¹)	75.8
Fluorescent Tag:	Coumarin
Cell Labeling:	Universal cell labeling reagent - No transfection agent required
Functionality	Amine
Visualization:	Florescence, Prussian blue, MRI, TEM
Storage Buffer:	dH ₂ O
Autoclaved:	No
Long-Term Storage:	4 °C - Do not freeze!
Similar Products:	None
Expiration Date:	None
